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<https://www.nature.com/articles/s41598-019-46875-5>

<https://swiftbiosci.com/publications/#methyl>

## Systematic evaluation of library preparation methods and sequencing platforms for high-throughput whole genome bisulfite sequencing

Li Zhou<sup>1</sup>, Hong Kiat Ng<sup>1</sup>, Daniela I. Drautz-Moses<sup>2</sup>, Stephan C. Schuster<sup>2</sup>, Stephan Beck<sup>3</sup>, Changhoon Kim<sup>4</sup>, John Campbell Chambers<sup>1,5,6,7</sup> & Marie Loh<sup>1,5,8</sup>

Swift  
BIOSCIENCES™

# THE GOLD STANDARD.

Swift's Accel-NGS® Methyl-Seq: Powered by Patented Adaptase® Technology!

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

\* CpX = CpG + CpH

Sample	Total Aligned	Genome Coverage	Duplicate Reads	Estimated Library Size	Relative Library Size	%CPX Missing	%CPX Covered >10x	
100 ng	<b>Swift</b>	<b>89.6%</b>	<b>22X</b>	<b>1.9%</b>	<b>713,928,657</b>	<b>1</b>	<b>0.6</b>	<b>92.2</b>
	Traditional	80.2%	21X	2.7%	603,864,885	0.85	0.6	88.1
	Random Primer	71.4%	16X	22.1%	47,841,839	0.07	7.7	39.4
10 ng	<b>Swift</b>	<b>87.8%</b>	<b>22X</b>	<b>2.7%</b>	<b>405,907,312</b>	<b>1</b>	<b>0.6</b>	<b>90.4</b>
	Traditional	76.8%	19X	11.9%	70,015,055	0.17	0.6	83.9
	Random Primer	71.9%	16X	22.2%	45,166,496	0.11	5.2	45.2
1 ng	<b>Swift</b>	<b>83.3%</b>	<b>18X</b>	<b>18.2%</b>	<b>38,295,868</b>	<b>1</b>	<b>0.6</b>	<b>77.1</b>
	Traditional	80.7%	10X	62.3%	6,215,803	0.16	2.0	17.0
	Random Primer	73.4%	12X	46.1%	11,863,382	0.31	6.6	31.3

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米国Swift Biosciences, Inc.  
<https://swiftbiosci.com/>

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## [Circulating tumor DNA profile recognizes transformation to castration-resistant neuroendocrine prostate cancer](#)

Loss of androgen receptor (AR) signaling dependence occurs in approximately 15%–20% of advanced treatment-resistant prostate cancers, and this may manifest clinically as transformation from a prostate adenocarcinoma histology to a castration-resistant neuroendocrine prostate cancer (CRPC-NE). (First published February 24, 2020)

# nature

## [SPEN integrates transcriptional and epigenetic control of X-inactivation](#)

*Xist* represents a paradigm for the function of long non-coding RNA in epigenetic regulation, although how it mediates X-chromosome inactivation (XCI) remains largely unexplained. (Published: 05 February 2020)

## SCIENTIFIC REPORTS

## [Ancient DNA and high-resolution chronometry reveal a long-term human role in the historical diversity and biogeography of the Bahamian hutia](#)

Quaternary paleontological and archaeological evidence often is crucial for uncovering the historical mechanisms shaping modern diversity and distributions. (Published: 28 January 2020)

bioRxiv  
THE PREPRINT SERVER FOR BIOLOGY

## [A versatile method for circulating cell-free DNA methylome profiling by reduced representation bisulfite sequencing](#)

The methylation profile of circulating cell-free DNA (cfDNA) in blood can be exploited to detect and diagnose cancer and other tissue pathologies and is therefore of great diagnostic interest. There is an urgent need for a cost-effective genome-wide methylation profiling method that is simple, robust and automatable and that works on highly fragmented cfDNA. (Posted June 11, 2019)

OXFORD  
ACADEMIC

## [Human mitochondrial DNA is extensively methylated in a non-CpG context](#)

Mitochondrial dysfunction plays critical roles in cancer development and related therapeutic response; however, exact molecular mechanisms remain unclear. (Published: 06 September 2019)

ScienceDirect

## [Epigenetic disturbances in obesity and diabetes: Epidemiological and functional insights](#)

Obesity and type 2 diabetes (T2D) are major public health issues worldwide, and put a significant burden on the healthcare system. Genetic variants, along with traditional risk factors such as diet and physical activity, could account for up to approximately a quarter of disease risk. (Available online 6 September 2019)