



Date: 2014/09/30

Information about Minor Change for CAX-371

Dear Customer:

We are pleased to inform that the following specifications had been changed for CAX-371.

Eco Mode

(Before) The chamber can cool down to the set temperature within approximately 3 minutes after release from ECO mode.

(After) The chamber can cool down to the set temperature within approximately 3 and half minutes after release from ECO mode.

Refrigerant

(Before) HFC 134a (280g)

(After) HFC R404A (320g)

Power Consumption (Heat output)

(Before) 1040W (900kcal / h)

(After) 930W (800kcal / h)

Net Weight

(Before) 118kg (120V) 123kg (220 / 230 / 240V)

(After) 115kg (120V) 120kg (220 / 230 / 240V)

Effected from September, 2014

If you have question, please contact for our staff.

TOMY DIGITAL BIOLOGY CO., LTD.

3-14-17 Tagara Nerima-kuTokyo 179-0073 Japan

phone: +81-3-5971-8160 fax: +81-3-3970-6036

<http://www.digital-biology.co.jp>



Hybrid Refrigerated Centrifuge
CAX-371