



Micro one

TOMY

Mini Personal Centrifuge

◀ Microtubes Rotor and Tube-strips Rotor Exchangeable. ▶



Spins 8xMicrotubes 1.5-2ml and 2x8-well Tube-strips.

The rotors can be installed with a single motion without tools.

This compact personal centrifuge will fit comfortably in the workspace.



User-friendly advanced design in attractive colors.

This made-in-Japan equipment promises you long lasting with highly reliable quality.

 **Digital Biology®**

Model	 microONE	
Rotor	 microtubes Rotor MA-1	 Tube-strips Rotor PA-1
Speed	6,000rpm (±20%)	6,200rpm (±20%)
Force	1,290 ~ 2,910G	Max.1,200 ~ 2,690G Min.810 ~ 1,820G
Maximum Volume	1.5(2.0)ml×8tubes	0.2ml×8tubes Connected ×2lines
Power Required	AC120V,1A,50/60Hz AC230V,1A,50/60Hz(CE)	
Fuse	250V,0.5A	
Dimensions	155W×175D×125Hmm	
Material	Lid	PC(UL94V-0)
	Body	ABS(UL94V-0)
	Actuator	ABS(UL94V-0)
	Rubber	NBR
	Rotor	ABS(UL94V-0)
Speed Control	Automatic / Constant Speed	
Drive System	Direct Drive / Vibration Insulated	
Brake	Dynamic Braking	
Acceleration Time	Within 5 seconds up to 90% of rated speed	
Deceleration Time	Within 15 seconds (Lid Close) Within 3 seconds (Lid Open)	

■ Accessories

Tube Adapter			
	For 0.5ml tube (Set of 8)		For 0.2ml tube (Set of 8)

Sales Office

TOMY DIGITAL BIOLOGY CO., LTD.

3-14-17 Tagara, Nerima-ku, Tokyo 179-0073, Japan
 e-mail: info@digital-biology.co.jp
 URL: http://www.digital-biology.co.jp
 phone: +81-3-5971-8160 fax: +81-3-3970-6036

Manufacturer

TOMY KOGYO CO., LTD.

3-14-17 Tagara, Nerima-ku, Tokyo 179-0073, Japan

TOMY TECH U.S.A., INC.

40479 Encyclopedia Circle, Fremont, California 94538, U.S.A.
 e-mail: info@tomytech.com
 URL: http://www.tomytech.com
 phone: 510-440-1976 fax: 510-440-1975
 Toll-Free US & Canada: 800-545-TOMY

TOMY SEIKO CO., LTD

All TOMY products have a limited one-year warranty.
 Specifications are subject to change according to product advancement.
 Tomy and Digital Biology is registered trademark of Tomy Seiko Co., Ltd.
 And Tomy Digital Biology Co., Ltd. Copyright 2008,
 Tomy Seiko and its subsidiaries. Printed in Japan.