

December 1, 2024

Notice of specification and model changes for Autoclave FLS-1000 Basket and Bucket

Dear Valued Customer:

Thank you for continued support for us.

This letter is to notify of specification changes for the following current products. The model will be changed according to the specification change.

Product	Current Model	Size (mm)	➔	New Model	Size (mm)	Autoclave
Basket	ACA-1000B	Φ415*233	➔	ACA-1100B	Φ415*236	FLS-1000
Bucket	ACA-1000K	Φ415*239	➔	ACA-1000KS	Φ415*239	

Basket

The specification changes are in the shape of the bottom stand and handle base part. When two baskets are stacked on top of each other, there is no space between the baskets and the total height is reduced.

Bucket

The specification change is in the handle base part. There is no change in size (inner dimension and height).

The current models will be discontinued, and the new models will be replaced as soon as the stock runs out. The basket packed with the Autoclave main unit as a standard accessory will also be replaced.

Contact:

TOMY DIGITAL BIOLOGY CO., LTD.

3-14-17, Tagara Nerima-ku, Tokyo 179-0073, Japan

Phone: +81-3-5971-8160

Email: marketing_mp@digital-biology.co.jp

<https://www.digital-biology.co.jp/manufactured/>