

October 1, 2025

Notice of model change

Continuous Flow Rotors

for High Speed Refrigerated Centrifuge Suprema25

Dear Valued Customer:

Thank you for your continued support of our products.

This letter is to notify the model change of the **Continuous Flow Rotors** for the High Speed Refrigerated Centrifuge Suprema25.

The previous models **NC-1**, **NC-2** and **NC-4** have been discontinued, and the new models **NCA-1200** and **NCC-1201** have been released.

<u>Previous Model</u>		<u>New Model</u>
Sealed aseptic type NC-2	➡	Sealed aseptic type NCA-1200
Collector type NC-1	➡	Collector type NCC-1201
Sealed aseptic type NC-4	➡	Discontinued

Suprema25 with serial numbers prior to **60389005** are not compatible with new models of continuous flow rotors. A new unit of Suprema25 is necessary to be obtained for use with the new continuous flow rotors.

For more product information, check below.

■New Continuous Flow Rotors: NCA-1200, NCC-1201

https://www.digital-biology.co.jp/manufactured/products/centrifuges/suprema/?abroad_anchor=sp_CFR

Contact:

TOMY DIGITAL BIOLOGY CO., LTD.

3-14-17, Tagara Nerima-ku, Tokyo 179-0073, Japan

Phone: +81-3-5971-8160

Email: marketing_mp@digital-biology.co.jp

<https://www.digital-biology.co.jp/manufactured/>