


ST-ZERO

Sterilizing Container for Wastes




If you have a hard time handling wastes in a biosafety cabinet, use our sterilizing container, ST-ZERO to store wastes inside the cabinet, and to allow sterilization with vertical autoclaves.



Key Features

- An easy to install design for a biosafety cabinet.
- Accommodating 25 X 25ml disposable pipettes.
- Stands provided to curb the effects on airflow. (See the installation examples on the right.)
- And leak proof while being sterilized in vertical autoclaves, as well as storing wastes within the cabinet.



Specifications

Material Stainless steel (SUS304)
 Dimensions(mm) . . . 106W,350D,80H(Protruding portion 380D)
 Weight 1.2kg
 Two stands included

Sales Office
TOMY DIGITAL BIOLOGY CO., LTD.
 3-14-17 Tagara, Nerima-ku, Tokyo 179-0073, Japan
 e-mail : info@digital-biology.co.jp
 URL : http://www.digital-biology.co.jp
 phone : +81-3-5971-8160 fax : +81-3-3970-6036

Manufacturer
TOMY KOGYO CO., LTD.
 3-14-17 Tagara, Nerima-ku, Tokyo 179-0073, Japan

TOMY SEIKO CO.,LTD

All TOMY products have a limited one-year warranty. Specifications are subject to change according to product advancement. Tomy and Digital Biology is registered trademark of Tomy Seiko Co., Ltd. And Tomy Digital Biology Co., Ltd. Copyright 2007, Tomy Seiko and its subsidiaries. Printed in Japan.