

■ SPECIFICATIONS

Single Tube Spin-able Centrifuge

Model	One Spin	
Max. Capacity	2.0 ml × 4 microtubes	
Tubes	2.0 ml × 4 microtubes 1.5 ml × 4 microtubes	
Speed	4,000rpm (±20%)	
RCF	450 ~1,000G (2.0 ml) 420 ~950G (1.5 ml)	
Safety Devices	Lid opening interlocking brake, Overcurrent fuse, Motor overcurrent detector (Motor lock, Unbalance)	
Compliance of Standard	EU Directives	
Dimensions	105W×130D×96H mm	
Net Weight	380 g	
Power Requirements	Dedicated AC adapter Input : 1-Phase AC100-240V 50/60Hz 0.6A	
Accessories Included	•AC adapter (6A-181WP05, 5 types of plugs) × 1set •Operator's manual × 1 copy	



- STANDARD ACCESSORY
- AC adapter
- 5 type power plug adapters

Cordless Multi-purpose Centrifuge

Model	Multi Spin		
Rotor	MA-1	PA-1	
Max. Capacity	2.0ml × 8 microtubes	0.2 ml × 8 PCR tubes (2 strips)	
Tubes	2.0 ml × 8 microtubes 1.5 ml × 8 microtubes	0.2 ml × 8 PCR tubes (2 strips)	
Speed	6,000rpm (±20%)	6,200rpm (±20%)	
RCF	1,310~2,960G (2.0 ml) 1,260~2,840G (1.5 ml)	Max. 1,200 ~2,690G Min. 810 ~1,820G	
Timer	6 minutes (Fixed)		
Safety Devices	Lid opening interlocking brake, Overcurrent fuse, Motor overcurrent detector (Motor lock, Unbalance)		
Compliance of Standard	EU Directives		
Dimensions	150W×163D×113H mm		
Net Weight	600 g (excluding rotors and batteries)		
Power Requirements	Dedicated AC adapter Input: 1-Phase AC100-240V 50/60Hz 0.6A or 4×AA size alkaline battery		
Accessories Included	AC adapter (6A-181WP05, 5 types of plugs) × 1 set Microtubes rotor MA-1 × 1pc. Tube-strips rotor PA-1 × 1pc. Operator's manual × 1copy		

Sales Office:

TOMY DIGITAL BIOLOGY CO., LTD.

3-14-17 Tagara, Nerima-ku, Tokyo 179-0073, Japan e-mail: info@digital-biology.co.jp
URL: http://www.digital-biology.co.jp/manufactured/phone: +81-3-5971-8160 fax: +81-3-3970-6036

TOMY SEIKO CO., LTD.

Manufacturer:

TOMY KOGYO CO., LTD.

3-14-17 Tagara, Nerima-ku, Tokyo 179-0073, Japan

All TOMY products have a limited one-year warranty.

Specifications are subject to change according to product advancement.


Tomy and Digital Biology is registered trademark of Tomy Seiko Co.,

Ltd. and Tomy Digital Biology Co., Ltd. Copyright 2016,

Tomy Seiko and its subsidiaries. Printed in Japan.







A single microtube is able to spin without balance load. (Patent pending)



Lid close to spin,
Release hand to stop

①Load tube ②Close lid ③Release (Start) hand(Sto







Easy tube unloading



Fits to universal power line


5 power plug adapters are included to fit to universal power line (AC100-240V). Comply with CE standard.




Spin-column or Screw capped tube available

Compact Design

Small foot print to save your bench space.







Motor over current & Imbalance detection system

6 minutes Timer equipped



• Do not use and leave centrifuge in extremely high/low temperature and in high humidity environment because it may cause serious damage to centrifuge.