

Expected Lifetime of Rotors and Buckets in Years or Cycles

Rotors and Buckets	Lifetime in years ^{*1}	Lifetime in cycles ^{*2}
NA-1, NA-3HS, NA-4HS, NA-6, NA-7, NA-8, NA-9, NA-11, NA-12, NA-16, NA-18, NA-20, NA-22, NA-22HS, NA-23	10	–
NA-400	10	–
NA-610	7	–
NC-1, NC-2, NC-4	10	–
NS-1, A29-95	10	–
TMP-1, TMP-11, TMP-21	10	50,000
TMP-2, TMP-24	10	50,000
TMA-22, TMA-26, TMA-27, TMA-29, TMA-30, TMA-32, TMA-25BH, TMA-100, TMA-200, TMA-300	10	–
CA-1, CA-10	10	50,000
CA-4HS, CA-5, CA-6, CA-8, CA-12, CA-14HS, CA-15, CA-16, CA-17, CA-100, CA-200, CA-300	10	–
CS-1, TMS-21, M0404-09, M0406-05, M0415-04	10	–
TS-4N, TS-4C, S4096-02	10	–
TS-6C, 6150-01, 6015-02	10	–
TS-33N, TS-33C, B433	10	–
TS-7N, TS-7C, 7050-02, 7150-01, 7050-01, 7115-08, 7015-08	10	–
7215-06, 7015-06, B407, SC-2, 7M5015-01	10	–
TS-38N, TS-38C	10	–
B438	10	40,000
B438-96	10	40,000
B438-29, B438-1507BH, B438-5002BH	10	–
TS-41N, TS-41C	10	–
B241	7	20,000
TS-36N	10	–
B436	10	53,000
TS-37C	10	–
B437	10	53,000
TS-12C, B212-1532, B212-0650	10	–

* Contact us for other rotors and buckets than the above listed.

*1 Period of time after dispatch from factory

*2 Number of times when used on the centrifuge independently