





Bottle


<ACA-700B>




500mL x 11pcs




1,000mL x 8pcs



2,000mL x 4pcs




5,000mL x 1pc




10,000mL x 1pc


<ACA-315B>




500mL x 8pcs




1,000mL x 5pcs



2,000mL x 3pcs



5,000mL x 1pc




10,000mL x 1pc




Flask


<ACA-700B>




100mL x 20pcs




300mL x 12pcs



500mL x 7pcs




1,000mL x 4pcs




5,000mL x 1pc


<ACA-315B>




100mL x 14pcs




300mL x 7pcs




500mL x 5pcs



1,000mL x 3pcs




5,000mL x 1pc



Animal Cage

<ACA-700B>




Internal ϕ 345mm

Internal 181mm

External 184mm

Loading Size




243x243mm

269x173mm

329x87mm

<ACA-315B>



Internal ϕ 300mm

Internal 181mm

External 190mm

Loading Size

210x210mm

258x151mm

288x76mm