


Rotors & Accessories Selection Guide

Continuous Flow Rotors for Suprema 25

1000ml




NC-1 [Click Here](#)
collector-type
1,000ml
17,000rpm : 31,350G



NC-2 [Click Here](#)
Germ-Free type
1,000ml
14,000rpm : 21,260G

NC-1

Applicable model	Suprema25
Max. speed	17,000 rpm
Max. RCF	31,350 G
Max. radius	97 mm
Min. radius	26 mm
Rotor capacity	1,000ml
Precipitation Capacity	500ml
Standard flow velocity	300ml/min



NC-2

Applicable model	Suprema25
Max. speed	14,000 rpm
Max. RCF	21,260 G
Max. radius	97 mm
Min. radius	26 mm
Rotor capacity	1,000ml
Precipitation Capacity	500ml
Standard flow velocity	300ml/min

