

# Rotors & Accessories Selection Guide

## Continuous Flow Rotors for Suprema 25

1,200ml



**NCA-1200**

[Click Here](#)

Sealed aseptic type  
Rotor capacity: 1,200ml  
Sedimentation capacity: 1,000ml  
17,000rpm 31,350G




**NCC-1201**

[Click Here](#)

Collector type  
Rotor capacity: 1,200ml  
Sedimentation capacity: 1,000ml  
17,000rpm 31,350G

## NCA-1200

Applicable model	Suprema25
Max. speed	17,000 rpm
Max. RCF	31,350 G
Max. radius	97 mm
Min. radius	23 mm
Rotor capacity	1,200ml
Sedimentation capacity	1,000ml
Standard flow velocity	300ml/min



## NCC-1201

Applicable model	Suprema25
Max. speed	17,000 rpm
Max. RCF	31,350 G
Max. radius	97 mm
Min. radius	23 mm
Rotor capacity	1,200ml
Sedimentation capacity	1,000ml
Standard flow velocity	300ml/min

